

## A FEW ADDED SUGGESTIONS



Do not twist your head as though to assist the ear, for the sound will strike the outside of the phone.



Face the speaker directly, so that the sound may enter the opening of the Super-Ear with full force.



The Super-Ear may be kept handy for instant use by hooking the headband around the collar as shown above.

For best results, start close up at first, and gradually increase the distance until you find your limit, according to your degree of deafness.

If both ears are equally impaired, an extra Super-Ear for the other ear will make an improvement in your results. This will cost \$10.00 extra including a double headband, and will also be subject to our trial guarantee.

The Super-Ear may be carried in the pocket or handbag, ready to be applied at once with or without the earloop.

It is usually some other member of the family who realizes FIRST what a help the Super-Ear is, as the speaker will not be required to strain his voice.

If the Super-Ear fits loosely on the ear when worn with the earloop, it may be further supported by a woman's hair, firmly pinned around it.

The earloop should press on top and bottom of the ear, and not in the back. You must practice how to adjust it. If the loop is not convenient in your case, you may use the headband.

**IMPORTANT:** Select the correct size ear tip to fit the orifice of your ear; large, small or medium. The ear tip is inserted by a push and a slight turn.

Bend the headband so that the Super-Ear and headband will fit snugly against the head.

Manufactured by

**American Earphone Co. Inc., N. Y.**