

Note: This card shows all accessories available. To order items that are not included as part of the standard kit, please contact your implant clinic or Cochlear. Images shown are not actual size.

1 ESPrIt™ 3G Speech Processor

Available in 4 colours:
Beige, Brown, Black, Silver.



2 ESPrIt 3G Coil and Cable

Available in 4 colours:
Beige, Brown, Black, Grey.
Available in 5 lengths:
6 cm, 8 cm, 10 cm,
20 cm, 28 cm.



3 Titanium Earhook

Can be bent to fit ear.



Reorder No. Z60118

4 Soft Earhook Cover

Can be worn over the titanium earhook for increased comfort.



Reorder No. Z60119

5 Earmould Compatible Earhook

For use with an earmould.



Reorder No. Z60120

6 Socket Cover

Attach to socket at base of processor when you are not using accessories.



Reorder No. Z60154

7 Drying Capsule

Put your ESPrIt 3G in its case every night with the drying capsule. The drying capsule will change colour when it needs replacing.



Reorder No. N94527

8 Cleaning Brush

Clean the battery contacts with the cleaning brush.



Reorder No. Z77029

9 Microphone Puffer

Use the microphone puffer to blow away moisture that may accumulate in the microphone ports.



Reorder No. Z60210

10 Batteries

We recommend Toshiba Custom High Power Zinc Air batteries which are available through Cochlear.



Reorder No. B13025

11 Screwdriver and Spare Screws

Use the screwdriver to remove the coil cover and change colours.



Reorder No. Z60209

12 Signal Check

Check if your ESPrIt 3G cable/coil is working.



Reorder No. Z22502

13 Accessories Adaptor

Attach to the base of the ESPrIt 3G to use accessories.



Reorder No. Z60121

14 Magnets and Spacers

A combination of magnets and spacers hold the transmitting coil over the implant.



Reorder No. Z60129 | Reorder No. Z60155

15 Accessory Adaptor Cable

Use to connect the TV/Hi-Fi cable and Personal Audio cable to the ESPril 3G.



Reorder No. Z77081

16 Coloured Battery and Coil Covers


12 vibrant colours that can be easily changed to suit your mood.



Reorder No. Z60073 Reorder No. Z60130

17 TV/Hi-Fi Cable

Connect a TV, stereo or personal computer to your ESPril 3G.



Reorder No. Z43135

18 Personal Audio Cable

Connect battery powered equipment, such as a portable CD player to your ESPril 3G.



Reorder No. Z43148

19 Lapel Microphone

Use to improve communication in noisy environments.



Reorder No. Z77127

20 Monitor Earphones

A hearing person can use these to check that sound is being received by the ESPril 3G microphone.



Reorder No. Z66133

Cochlear Ltd (ABN 96 002 618 073) 14 Mars Road, Lane Cove NSW 2066, Australia Tel: 61 2 9428 6555 Fax: 61 2 9428 6352
Cochlear Corporation 400 Inverness Drive South, Suite 400, Englewood CO 80112, USA Tel: 1 303 790 9010 Fax: 1 303 792 9025
Cochlear Europe Ltd 22-24 Worpole Road, Wimbledon, London SW19 4DD, United Kingdom Tel: 44 20 8879 4900 Fax: 44 20 8946 9066
Cochlear AG Margarethenstrasse 47, CH-4053 Basel, Switzerland Tel: 41 61 205 0404 Fax: 41 61 205 0405
Cochlear GmbH Karl-Wiechert-Allee 76A, D-30625 Hannover, Germany Tel: 49 511 542 770 Fax: 49 511 542 7770
Nihon Cochlear Co Ltd Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 0245
Cochlear (HK) Ltd 21F Shun Ho Tower, 24-30 Ice House Street, Central, Hong Kong Tel: 852 2530 5773 Fax: 852 2530 5183

ESPril is a trademark of Cochlear Limited.
PRINTED IN AUSTRALIA

The Cochlear logo, featuring three stylized, overlapping curved lines above the word "Cochlear" in a bold, sans-serif font.

Cochlear

N94B20F ISS2 JAN2002