


ESPril™ 3G Jacket




The jacket protects your ESPril™ 3G speech processor from moisture and humidity. You may find it useful if you are involved in active sports, live in a humid climate or perspire heavily.

The jacket should be worn every day. Remove it at the end of each day, prior to storing your ESPril 3G with an activated drying capsule in the ESPril 3G case.

With titanium and earmould compatible earhooks*



With the soft earhook cover on the earhook



Cochlear Ltd (ABN 96 002 818 073) 14 Mars Road, Lane Cove NSW 2066, Australia Tel: 61 2 9428 6555 Fax: 61 2 9438 6332
Cochlear Americas 400 Inverness Parkway, Suite 400, Englewood CO 80112, USA Tel: 1 303 790 9010 Fax: 1 303 792 9025
Cochlear Europe Ltd 22-24 Worpole Road, Wimborne, Dorset SW19 4DD, United Kingdom Tel: 44 20 8879 4900 Fax: 44 20 8946 9066
Cochlear AG Pflanzstrasse 47, CH - 4053 Basel, Switzerland Tel: 41 61 205 0404 Fax: 41 61 205 0405
Cochlear GmbH Karl-Wilhelm-Allee 76A, D-30425 Hannover, Germany Tel: 49 51 1 542 770 Fax: 49 51 1 542 7770
Nihon Cochlear Co Ltd Ochanomizu-Hotomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 0245
Cochlear (HK) Ltd 21F Shun Ho Tower, 24-30 Ice House Street, Central, Hong Kong Tel: 852 2530 5773 Fax: 852 2530 5183

www.cochlear.com



Nucleus is a registered trademark of Cochlear Limited.
ESPrIt is a trademark of Cochlear Limited.

PRINTED IN AUSTRALIA
N95540F ISS1